



The benefits of custom orthoses by ORTEMA

We manufacture custom orthoses for your patients. The carbon fiber construction ensures maximum comfort with minimal weight and maximum stability.

Our orthoses for the knee, shoulder, elbow, hand and face undergo strict quality controls and are the result of continuous research and development. This keeps us at the forefront of medical innovation and guarantees our customers the optimal solution.

Your benefits at a glance:

Precision Fit

Each orthosis is made according to individual requirements to ensure perfect fit and maximum comfort.

Functional Support

Our orthoses provide targeted relief and functional support, effectively helping in rehabilitation for injuries and providing preventive protection.

High-Quality Materials

We use only materials that meet the highest quality standards and are hypoallergenic to avoid skin irritation.

Versatile Applications

Whether it's the knee, wrist, elbow, shoulder, or head - ORTEMA offers the right solution for every need.



Supply Your Patients with the Proven ORTEMA Technique

ORTEMA GmbH Headquarter

Kurt-Lindemann-Weg 10
71706 Markgröningen
+49 7145 - 91 53 800
+49 7145 - 91 53 980
info@ortema.de

ortema.de



Version 05/2024



Orthopaedic Technology & Medical Supply Store

Supply your Patients with the ORTEMA Technique

Our Service Production for Medical Supply Stores and Orthopedic Stores



Discover the proven quality of ORTEMA for your patients. With our service manufacturing, we offer orthopedic technicians and medical supply stores the opportunity to have custom-made orthoses manufactured according to individual specifications, crafted to the highest standards of quality.

And it's as simple as that: You send us the scan data of your patient along with the documentation, and we will deliver a ready-to-try orthosis within ten working days, plus postal delivery time.

Benefit from the advantages of the ORTEMA Service Production:

- ▶ Time and cost efficiency
- ▶ Bridging staff shortages
- ▶ Developing new business fields
- ▶ Benefiting from our expertise

Have we sparked your interest?

Please request the ORTEMA Service Production documents by phone, fax or email:

Phone: +49 7145 - 91 53 800

Fax: +49 7145 - 91 53 980

email: servicefertigung@ortema.de

GET YOUR INDIVIDUAL ORTEMA ORTHOSIS IN JUST 5 EASY STEPS

(using the example of a K-COM knee brace)



1. Order

If you have a patient in need of a custom K-COM knee orthosis, please contact us to initiate the Service Production order.

2. 3D-Scan

You create a 3D scan of the patient's leg according to our specifications and send the data to us via email or upload.

3. Manufacturing

We create the 3D model and manufacture the K-COM knee orthosis according to the exact indication.

4. Adjustment

The trial-ready K-COM knee orthosis is then sent to your medical supply store or orthopedic store.

5. Inspection

Following the fitting, the treating physician can conduct a functional test of the orthosis to ensure high-quality care.

The production time is within 10 working days after receiving the scan, plus shipping time.